Static GHG Chambers



PUSH ON LID 250mm diameter



1/4 TURN LID 250mm diameter



PVC BASE RING + CHAMBER 250mm diameter Ring stands 30mm above soil



CHAMBER + EXTENSION + BASE 500x500mm. 250mm high chamber with fan, vent and sampling port. 500mm high extension. Stainless steel base.



VACUUM MANIFOLD Extraction of multiple vials for sample collection when connected to vacuum pump

Gibson Engineering specialises in research and sampling equipment. We can develop and customise all products to meet your needs. We specialise in design and manufacture of products from chambers to soil sampling equipment, laboratory equipment and anything else you may require.





sales@gibsonengineering.com.au 0405 041 654