

Features

PLC controlled operating system consisting of variable speed drive, countdown timer, time start/stop set points sealed in force vented compartment

Safety switch in lid to prevent machine operation when open

Emergency Stop button on machine front (not shown)

Smooth maintenance free drive train consisting of modern drive train components including quality gearbox and motor sealed in separate force vented compartment for clean trouble free operation.

Polished stainless steel external machine shroud for maintenance free use with spillages of KCl etc

Fully sealed lid – machine joint to decrease noise in laboratory from rotating elements in jars

50 Sample capacity in 100mL square jars